

## SEN2 Polyclonal Antibody

Cat No: HR1AP11508

For research use only

### Overview

Product Name	SEN2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	51kDa
GeneID?Human?	TSEN2 SEN2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	tRNA splicing endonuclease subunit 2(TSEN2) Homo sapiens This gene encodes one of the subunits of the tRNA splicing endonuclease. This endonuclease catalyzes the first step in RNA splicing which is the removal of introns. Mutations in this gene have been associated with pontocerebellar hypoplasia type 2. A pseudogene has been identified on chromosome 4. Multiple transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Feb 2009],