

## EI2BB Polyclonal Antibody

Cat No: HR1AP11503

For research use only

### Overview

Product Name	EI2BB Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	38kDa
GeneID?Human?	EIF2B2 EIF2BB
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	eukaryotic translation initiation factor 2B subunit beta(EIF2B2) Homo sapiens This gene encodes the beta subunit of eukaryotic initiation factor-2B (EIF2B). EIF2B is involved in protein synthesis and exchanges GDP and GTP for its activation and deactivation. [provided by RefSeq, Aug 2011],