

# SCRT1 Polyclonal Antibody

Cat No: HR1AP11498

For research use only

## Overview

Product Name	SCRT1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	38kDa
GeneID?Human?	SCRT1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	scratch family transcriptional repressor 1(SCRT1) Homo sapiens This gene encodes a C2H2-type zinc finger transcriptional repressor that binds to E-box motifs. The encoded protein may promote neural differentiation and may be involved in cancers with neuroendocrine features. [provided by RefSeq, Jul 2013],