

TRRAP Polyclonal Antibody

Cat No: HR1AP11499

For research use only

Overview

Product Name	TRRAP Polyclonal Antibody
Source	Rabbit
Applications	IHC-p
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	424kDa
GeneID?Human?	TRRAP PAF400
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>transformation/transcription domain associated protein(TRRAP) Homo sapiens This gene encodes a large multidomain protein of the phosphoinositide 3-kinase-related kinases (PIKK) family. The encoded protein is a common component of many histone acetyltransferase (HAT) complexes and plays a role in transcription and DNA repair by recruiting HAT complexes to chromatin. Deregulation of this gene may play a role in several types of cancer including glioblastoma multiforme. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Sep 2011],</p>