

# TTF2 Polyclonal Antibody

Cat No: HR1AP11494

For research use only

## Overview

Product Name	TTF2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	127kDa
GeneID?Human?	TTF2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	transcription termination factor 2(TTF2) Homo sapiens This gene encodes a member of the SWI2/SNF2 family of proteins, which play a critical role in altering protein-DNA interactions. The encoded protein has been shown to have dsDNA-dependent ATPase activity and RNA polymerase II termination activity. This protein interacts with cell division cycle 5-like, associates with human splicing complexes, and plays a role in pre-mRNA splicing. [provided by RefSeq, Jul 2008].