

TAF1L Polyclonal Antibody

Cat No: HR1AP11489

For research use only

Overview

Product Name	TAF1L Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	200kDa
GeneID?Human?	TAF1L
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	TATA-box binding protein associated factor 1 like(TAF1L) Homo sapiens This locus is intronless, and apparently arose in the primate lineage from retrotransposition of the transcript from the multi-exon TAF1 locus on the X chromosome. The gene is expressed in male germ cells, and the product has been shown to function interchangeably with the TAF1 product. [provided by RefSeq, Aug 2015],