



# LBX1 Polyclonal Antibody

Cat No: HR1AP11487

For research use only

## Overview

Product Name	LBX1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	30kDa
GeneID?Human?	LBX1 LBX1H
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>ladybird homeobox 1(LBX1) Homo sapiens This gene and the orthologous mouse gene were found by their homology to the Drosophila lady bird early and late homeobox genes. In the mouse, this gene is a key regulator of muscle precursor cell migration and is required for the acquisition of dorsal identities of forelimb muscles. [provided by RefSeq, Jul 2008],</p>