

BTF3 Polyclonal Antibody

Cat No: HR1AP11485

For research use only

Overview

Product Name	BTF3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	22kDa
GeneID?Human?	BTF3 NACB OK/SW-cl.8
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	basic transcription factor 3(BTF3) Homo sapiens This gene encodes the basic transcription factor 3. This protein forms a stable complex with RNA polymerase IIB and is required for transcriptional initiation. Alternative splicing results in multiple transcript variants encoding different isoforms. This gene has multiple pseudogenes. [provided by RefSeq, Jul 2008],