

TCF25 Polyclonal Antibody

Cat No: HR1AP11484

For research use only

Overview

Product Name	TCF25 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	74kDa
GeneID?Human?	TCF25 KIAA1049 NULP1 FKSG26
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	transcription factor 25(TCF25) Homo sapiens TCF25 is a member of the basic helix-loop-helix (bHLH) family of transcription factors that are important in embryonic development (Steen and Lindholm, 2008 [PubMed 18068114]).[supplied by OMIM, Sep 2008],