

TCF20 Polyclonal Antibody

Cat No: HR1AP11483

For research use only

Overview

Product Name	TCF20 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	215kDa
GeneID?Human?	TCF20 KIAA0292 SPBP
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	transcription factor 20(TCF20) Homo sapiens This gene encodes a transcription factor that recognizes the platelet-derived growth factor-responsive element in the matrix metalloproteinase 3 promoter. The encoded protein is thought to be a transcriptional coactivator, enhancing the activity of transcription factors such as JUN and SP1. Mutations in this gene are associated with autism spectrum disorders. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2015],