

# TCF15 Polyclonal Antibody

Cat No: HR1AP11482

For research use only

## Overview

Product Name	TCF15 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	21kDa
GeneID?Human?	TCF15 BHLHA40 BHLHEC2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	transcription factor 15 (basic helix-loop-helix)(TCF15) Homo sapiens The protein encoded by this gene is found in the nucleus and may be involved in the early transcriptional regulation of patterning of the mesoderm. The encoded basic helix-loop-helix protein requires dimerization with another basic helix-loop-helix protein for efficient DNA binding. [provided by RefSeq, Jul 2008],