

TCAL5 Polyclonal Antibody

Cat No: HR1AP11479

For research use only

Overview

Product Name	TCAL5 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	22kDa
GeneID?Human?	TCEAL5
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	transcription elongation factor A like 5(TCEAL5) Homo sapiens This gene, which is located on the X chromosome, encodes a protein which contains a BEX (brain expressed X-linked like family) domain. This domain is found in proteins encoded by the TCEAL elongation factor (transcription elongation factor A (SII)-like) gene family also located on the X chromosome. The coding region for this gene is located entirely in the terminal exon. [provided by RefSeq, Sep 2011],