

VGLL2 Polyclonal Antibody

Cat No: HR1AP11476

For research use only

Overview

Product Name	VGLL2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	34kDa
GeneID?Human?	VGLL2 VITO1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>vestigial like family member 2(VGLL2) Homo sapiens This gene encodes a protein with a transcriptional enhancer factor 1 (TEF-1) interaction domain. The encoded protein may act as a co-factor of TEF-1 regulated gene expression during skeletal muscle development. Alternatively spliced transcript variants encoding multiple isoforms have been observed. [provided by RefSeq, Jul 2012],</p>