

## TPC6B Polyclonal Antibody

Cat No: HR1AP11475

For research use only

### Overview

Product Name	TPC6B Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	17kDa
GeneID?Human?	TRAPPC6B
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	trafficking protein particle complex 6B(TRAPPC6B) Homo sapiens TRAPPC6B is a component of TRAPP complexes, which are tethering complexes involved in vesicle transport (Kummel et al., 2005 [PubMed 16025134]).[supplied by OMIM, Mar 2008],