

# VGLL4 Polyclonal Antibody

Cat No: HR1AP11477

For research use only

## Overview

Product Name	VGLL4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	31kDa
GeneID?Human?	VGLL4 KIAA0121
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:May act as a specific coactivator for the mammalian TEFs.,similarity:Belongs to the vestigial family.,subunit:Interacts with TEFs.,