

TPPC4 Polyclonal Antibody

Cat No: HR1AP11474

For research use only

Overview

Product Name	TPPC4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	24kDa
GeneID?Human?	TRAPPC4 SBDN CGI-104 HSPC172 PTD009
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:May play a role in vesicular transport from endoplasmic reticulum to Golgi.,similarity:Belongs to the TRAPP small subunits family. TRAPPC4 subfamily.,subunit:Part of the multisubunit TRAPP (transport protein particle) complex. Interacts with SDC2.,