

TRAK1 Polyclonal Antibody

Cat No: HR1AP11472

For research use only

Overview

Product Name	TRAK1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	104kDa
GeneID?Human?	TRAK1 KIAA1042 OIP106
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	PTM:O-glycosylated.,similarity:Contains 1 HAP1 N-terminal domain.,subunit:Interacts with O-GlcNAc transferase. Interacts with RHOT1/Miro-1 and RHOT2/Miro-2.,tissue specificity:High expression in spinal cord and moderate expression in all other tissues and specific brain regions examined. Expressed in all cell lines examined.,