

TM1L1 Polyclonal Antibody

Cat No: HR1AP11471

For research use only

Overview

Product Name	TM1L1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	52kDa
GeneID?Human?	TOM1L1 SRCASM
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>function:Probable adapter protein involved in signaling pathways. Interacts with the SH2 and SH3 domains of various signaling proteins when it is phosphorylated. May promotes FYN activation, possibly by disrupting intramolecular SH3-dependent interactions.,PTM:Phosphorylated on tyrosines by FYN (By similarity). Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Belongs to the TOM1 family.,similarity:Contains 1 GAT domain.,similarity:Contains 1 VHS domain.,subunit:Interacts with FYN, GRB2 and PIK3R1 when phosphorylated.,</p>