

## TIPRL Polyclonal Antibody

Cat No: HR1AP11469

For research use only

### Overview

Product Name	TIPRL Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	29kDa
GeneID?Human?	TIPRL
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	TOR signaling pathway regulator(TIPRL) Homo sapiens TIPRL is an inhibitory regulator of protein phosphatase-2A (PP2A) (see PPP2CA; MIM 176915), PP4 (see PPP4C; MIM 602035), and PP6 (see PPP6C; MIM 612725) (McConnell et al., 2007 [PubMed 17384681]).[supplied by OMIM, Nov 2010],