

TENN Polyclonal Antibody

Cat No: HR1AP11461

For research use only

Overview

Product Name	TENN Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	142kDa
GeneID?Human?	TNN
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Involved in neurite outgrowth and cell migration in hippocampal explants.,similarity:Belongs to the tenascin family.,similarity:Contains 1 fibrinogen C-terminal domain.,similarity:Contains 3 EGF-like domains.,similarity:Contains 9 fibronectin type-III domains.,