

TPRX1 Polyclonal Antibody

Cat No: HR1AP11462

For research use only

Overview

Product Name	TPRX1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	45kDa
GeneID?Human?	TPRX1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	tetrapeptide repeat homeobox 1 (TPRX1) Homo sapiens Homeobox genes encode DNA-binding proteins, many of which are thought to be involved in early embryonic development. Homeobox genes encode a DNA-binding domain of 60 to 63 amino acids referred to as the homeodomain. This gene is a member of the TPRX homeobox gene family. [provided by RefSeq, Jul 2008],