

## TLX3 Polyclonal Antibody

Cat No: HR1AP11458

For research use only

### Overview

Product Name	TLX3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	32kDa
GeneID?Human?	TLX3 HOX11L2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	T-cell leukemia homeobox 3(TLX3) Homo sapiens The protein encoded by this gene is an orphan homeobox protein that encodes a DNA-binding nuclear transcription factor. A translocation [t(5;14)(q35;q32)] involving this gene is associated with T-cell acute lymphoblastic leukemia (T-ALL) in children and young adults. [provided by RefSeq, Nov 2015],