

TBX15 Polyclonal Antibody

Cat No: HR1AP11456

For research use only

Overview

Product Name	TBX15 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	66kDa
GeneID?Human?	TBX15 TBX14
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	T-box 15(TBX15) Homo sapiens This gene belongs to the T-box family of genes, which encode a phylogenetically conserved family of transcription factors that regulate a variety of developmental processes. All these genes contain a common T-box DNA-binding domain. Mutations in this gene are associated with Cousin syndrome.[provided by RefSeq, Oct 2009],