

TBD2A Polyclonal Antibody

Cat No: HR1AP11455

For research use only

Overview

| | |
|-----------------------|--|
| Product Name | TBD2A Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,ELISA |
| Species Reactivity | Human |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 102kDa |
| GeneID?Human? | TBC1D2 PARIS1 PP8997 TBC1D2A |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | |
| Background | function:May act as a GTPase-activating protein for Rab family protein(s).,similarity:Contains 1 PH domain.,similarity:Contains 1 Rab-GAP TBC domain.,tissue specificity:Expressed in a broad range of tissues. Differentially expressed in prostate normal and cancer cells., |