

STX6 Polyclonal Antibody

Cat No: HR1AP11453

For research use only

Overview

Product Name	STX6 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	28kDa
GeneID?Human?	STX6
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Involved in intracellular vesicle trafficking.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Belongs to the syntaxin family.,similarity:Contains 1 t-SNARE coiled-coil homology domain.,subunit:Binds EEA1. Interacts with VPS45A. Interacts with MARCH2 and MARCH3 (By similarity). Interacts with GOPC. Identified in a complex containing STX6, STX13, VAMP4 and VT11A.,