



## SNPH Polyclonal Antibody

Cat No: HR1AP11451

For research use only

### Overview

Product Name	SNPH Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	54kDa
GeneID?Human?	SNPH KIAA0374
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>syntaphilin(SNPH) Homo sapiens Syntaxin-1, synaptobrevin/VAMP, and SNAP25 interact to form the SNARE complex, which is required for synaptic vesicle docking and fusion. The protein encoded by this gene is membrane-associated and inhibits SNARE complex formation by binding free syntaxin-1. Expression of this gene appears to be brain-specific. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Dec 2015],</p>