

SYT15 Polyclonal Antibody

Cat No: HR1AP11444

For research use only

Overview

Product Name	SYT15 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	46kDa
GeneID?Human?	SYT15
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>synaptotagmin 15(SYT15) Homo sapiens This gene encodes a member of the Synaptotagmin (Syt) family of membrane trafficking proteins. Members of this family contain a transmembrane region and a C-terminal-type tandem C2 domain. Unlike related family members, the encoded protein may be involved in membrane trafficking in non-neuronal tissues. Two transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq, Jul 2008].</p>