

SV2A Polyclonal Antibody

Cat No: HR1AP11439

For research use only

Overview

Product Name	SV2A Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	81kDa
GeneID?Human?	SV2A KIAA0736 PSEC0174
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	synaptic vesicle glycoprotein 2A(SV2A) Homo sapiens The protein encoded by this gene is one of three related synaptic vesicle proteins. The encoded protein may interact with synaptotagmin to enhance low frequency neurotransmission in quiescent neurons. [provided by RefSeq, Jun 2016],