

SIKE1 Polyclonal Antibody

Cat No: HR1AP11434

For research use only

Overview

Product Name	SIKE1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	22kDa
GeneID?Human?	SIKE1 SIKE
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>suppressor of IKBKE 1(SIKE1) Homo sapiens SIKE interacts with IKK-epsilon (IKBKE; MIM 605048) and TBK1 (MIM 604834) and acts as a suppressor of TLR3 (MIM 603029) and virus-triggered interferon activation pathways (Huang et al., 2005 [PubMed 16281057]).[supplied by OMIM, Mar 2008],</p>