



Z280B Polyclonal Antibody

Cat No: HR1AP11433

For research use only

Overview

Product Name	Z280B Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	59kDa
GeneID?Human?	ZNF280B SUHW2 ZNF279 ZNF632
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>zinc finger protein 280B(ZNF280B) Homo sapiens The protein encoded by this gene is a transcription factor that upregulates expression of MDM2, which negatively regulates p53 expression. This gene is highly expressed in prostate cancer cells, which leads to a reduction in p53 levels and an increase in growth of the cancer cells. Several transcript variants have been found for this gene, but only one of them is protein-coding. [provided by RefSeq, Jan 2015],</p>