

SOCS5 Polyclonal Antibody

Cat No: HR1AP11431

For research use only

Overview

Product Name	SOCS5 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	58kDa
GeneID?Human?	SOCS5 CIS6 CISH5 CISH6 KIAA0671
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>suppressor of cytokine signaling 5(SOCS5) Homo sapiens The protein encoded by this gene contains a SH2 domain and a SOCS BOX domain. The protein thus belongs to the suppressor of cytokine signaling (SOCS) family, also known as STAT-induced STAT inhibitor (SSI) protein family. SOCS family members are known to be cytokine-inducible negative regulators of cytokine signaling. The specific function of this protein has not yet been determined. Two alternatively spliced transcript variants encoding an identical protein have been reported. [provided by RefSeq, Jul 2008].</p>