

SOCS4 Polyclonal Antibody

Cat No: HR1AP11430

For research use only

Overview

Product Name	SOCS4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	48kDa
GeneID?Human?	SOCS4 SOCS7
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>suppressor of cytokine signaling 4(SOCS4) Homo sapiens The protein encoded by this gene contains a SH2 domain and a SOCS BOX domain. The protein thus belongs to the suppressor of cytokine signaling (SOCS), also known as STAT-induced STAT inhibitor (SSI), protein family. SOCS family members are known to be cytokine-inducible negative regulators of cytokine signaling. Two alternatively spliced transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Jul 2008],</p>