



# SPEG Polyclonal Antibody

Cat No: HR1AP11428

For research use only

## Overview

Product Name	SPEG Polyclonal Antibody
Source	Rabbit
Applications	IHC-p
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	359kDa
GeneID?Human?	SPEG APEG1 KIAA1297
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	SPEG complex locus(SPEG) Homo sapiens This gene encodes a protein with similarity to members of the myosin light chain kinase family. This protein family is required for myocyte cytoskeletal development. Along with the desmin gene, expression of this gene may be controlled by the desmin locus control region. Mutations in this gene are associated with centronuclear myopathy 5. [provided by RefSeq, Jun 2016],