

STAR Polyclonal Antibody

Cat No: HR1AP11426

For research use only

Overview

Product Name	STAR Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	31kDa
GeneID?Human?	STAR STARD1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	steroidogenic acute regulatory protein(STAR) Homo sapiens The protein encoded by this gene plays a key role in the acute regulation of steroid hormone synthesis by enhancing the conversion of cholesterol into pregnenolone. This protein permits the cleavage of cholesterol into pregnenolone by mediating the transport of cholesterol from the outer mitochondrial membrane to the inner mitochondrial membrane. Mutations in this gene are a cause of congenital lipoid adrenal hyperplasia (CLAH), also called lipoid CAH. A pseudogene of this gene is located on chromosome 13. [provided by RefSeq, Jul 2008],