

# YSK4 Polyclonal Antibody

Cat No: HR1AP11424

For research use only

## Overview

Product Name	YSK4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	146kDa
GeneID?Human?	YSK4 RCK
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	catalytic activity:ATP + a protein = ADP + a phosphoprotein.,similarity:Belongs to the protein kinase superfamily. STE Ser/Thr protein kinase family. STE20 subfamily.,similarity:Contains 1 protein kinase domain.,