

## CLASR Polyclonal Antibody

Cat No: HR1AP11421

For research use only

### Overview

Product Name	CLASR Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	74kDa
GeneID?Human?	CLASRP SFRS16 SWAP2 UNQ2428/PRO4988
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	caution:It is uncertain whether Met-1 or Met-16 is the initiator.,function:Probably functions as an alternative splicing regulator. May regulate the mRNA splicing of genes such as CLK1. May act by regulating members of the CLK kinase family.,PTM:Phosphorylated in vitro by CLK4.,similarity:Belongs to the splicing factor SR family.,subunit:Probably interacts with CLK4.,