

## SREK1 Polyclonal Antibody

Cat No: HR1AP11420

For research use only

### Overview

Product Name	SREK1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	55kDa
GeneID?Human?	SREK1 SFRS12 SRRP86
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>splicing regulatory glutamic acid and lysine rich protein 1(SREK1) Homo sapiens This gene encodes a member of a family of serine/arginine-rich (SR) splicing proteins containing RNA recognition motif (RRM) domains. The encoded protein interacts with other SR proteins to modulate splice site selection. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Jul 2012].</p>