

SPF45 Polyclonal Antibody

Cat No: HR1AP11419

For research use only

Overview

Product Name	SPF45 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	44kDa
GeneID?Human?	RBM17 SPF45
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	RNA binding motif protein 17(RBM17) Homo sapiens This gene encodes an RNA binding protein. The encoded protein is part of the spliceosome complex and functions in the second catalytic step of mRNA splicing. Alternatively spliced transcript variants have been described. Related pseudogenes exist on chromosomes 9 and 15. [provided by RefSeq, Mar 2009],