

SUGP1 Polyclonal Antibody

Cat No: HR1AP11418

For research use only

Overview

Product Name	SUGP1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	70kDa
GeneID?Human?	SUGP1 SF4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	SURP and G-patch domain containing 1(SUGP1) Homo sapiens SF4 is a member of the SURP family of splicing factors.[supplied by OMIM, Sep 2003],