

SF3A1 Polyclonal Antibody

Cat No: HR1AP11416

For research use only

Overview

Product Name	SF3A1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	87kDa
GeneID?Human?	SF3A1 SAP114
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>splicing factor 3a subunit 1(SF3A1) Homo sapiens This gene encodes a subunit of the splicing factor 3a protein complex. The splicing factor 3a heterotrimer is a component of the mature U2 small nuclear ribonucleoprotein particle (snRNP). U2 small nuclear ribonucleoproteins play a critical role in spliceosome assembly and pre-mRNA splicing. [provided by RefSeq, Aug 2014],</p>