

SPAG5 Polyclonal Antibody

Cat No: HR1AP11415

For research use only

Overview

Product Name	SPAG5 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	131kDa
GeneID?Human?	SPAG5
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	sperm associated antigen 5(SPAG5) Homo sapiens This gene encodes a protein associated with the mitotic spindle apparatus. The encoded protein may be involved in the functional and dynamic regulation of mitotic spindles. [provided by RefSeq, Jul 2008],