

## SAS10 Polyclonal Antibody

Cat No: HR1AP11412

For research use only

### Overview

Product Name	SAS10 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	52kDa
GeneID?Human?	UTP3 CRLZ1 SAS10
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Essential for gene silencing: has a role in the structure of silenced chromatin. Plays a role in the developing brain.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Belongs to the SAS10 family.,