

S12A2 Polyclonal Antibody

Cat No: HR1AP11409

For research use only

Overview

Product Name	S12A2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	133kDa
GeneID?Human?	SLC12A2 NKCC1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>solute carrier family 12 member 2(SLC12A2) Homo sapiens The protein encoded by this gene mediates sodium and chloride transport and reabsorption. The encoded protein is a membrane protein and is important in maintaining proper ionic balance and cell volume. This protein is phosphorylated in response to DNA damage. Three transcript variants encoding two different isoforms have been found for this gene. [provided by RefSeq, Jan 2012],</p>