

SC6A7 Polyclonal Antibody

Cat No: HR1AP11408

For research use only

Overview

Product Name	SC6A7 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	69kDa
GeneID?Human?	SLC6A7 PROT
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>solute carrier family 6 member 7(SLC6A7) Homo sapiens This gene is a member of the gamma-aminobutyric acid (GABA) neurotransmitter gene family and encodes a high-affinity mammalian brain L-proline transporter protein. This transporter protein differs from other sodium-dependent plasma membrane carriers by its pharmacological specificity, kinetic properties, and ionic requirements. [provided by RefSeq, Jul 2008],</p>