

S20A1 Polyclonal Antibody

Cat No: HR1AP11407

For research use only

Overview

Product Name	S20A1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	74kDa
GeneID?Human?	SLC20A1 GLVR1 PIT1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>solute carrier family 20 member 1(SLC20A1) Homo sapiens The protein encoded by this gene is a sodium-phosphate symporter that absorbs phosphate from interstitial fluid for use in cellular functions such as metabolism, signal transduction, and nucleic acid and lipid synthesis. The encoded protein is also a retroviral receptor, causing human cells to be susceptible to infection by gibbon ape leukemia virus, simian sarcoma-associated virus, feline leukemia virus subgroup B, and 10A1 murine leukemia virus.[provided by RefSeq, Mar 2011],</p>