

NPT2A Polyclonal Antibody

Cat No: HR1AP11405

For research use only

Overview

| | |
|-----------------------|--|
| Product Name | NPT2A Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,ELISA |
| Species Reactivity | Human,Rat,Mouse |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 70kDa |
| GeneID?Human? | SLC34A1 NPT2 SLC17A2 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | |
| Background | solute carrier family 34 member 1(SLC34A1) Homo sapiens This gene encodes a member of the type II sodium-phosphate cotransporter family. Mutations in this gene are associated with hypophosphatemia nephrolithiasis/osteoporosis 1. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Nov 2009], |