

## S6A19 Polyclonal Antibody

Cat No: HR1AP11404

For research use only

### Overview

Product Name	S6A19 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	69kDa
GeneID?Human?	SLC6A19 B0AT1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	solute carrier family 6 member 19(SLC6A19) Homo sapiens This gene encodes a system B(0) transmembrane protein that actively transports most neutral amino acids across the apical membrane of epithelial cells. Mutations in this gene result in Hartnup disorder. [provided by RefSeq, Jul 2008],