

NAGT1 Polyclonal Antibody

Cat No: HR1AP11403

For research use only

Overview

Product Name	NAGT1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	56kDa
GeneID?Human?	NAGLT1 KIAA1919 HSPC100
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:May function as a sodium-dependent glucose transporter.,similarity:Belongs to the major facilitator superfamily.,