

SL9A5 Polyclonal Antibody

Cat No: HR1AP11397

For research use only

Overview

Product Name	SL9A5 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	H,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	98kDa
GeneID?Human?	SLC9A5 NHE5
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>caution:The number, localization and denomination of hydrophobic domains in the Na(+)/H(+) exchangers vary among authors.,function:Involved in pH regulation to eliminate acids generated by active metabolism or to counter adverse environmental conditions. Major proton extruding system driven by the inward sodium ion chemical gradient. Plays an important role in signal transduction.,PTM:Phosphorylated (Possible),similarity:Belongs to the monovalent cation:proton antiporter 1 (CPA1) transporter (TC 2.A.36) family.,tissue specificity:Expressed in brain, testis, spleen, and skeletal muscle.,</p>