

SL9A2 Polyclonal Antibody

Cat No: HR1AP11396

For research use only

Overview

Product Name	SL9A2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	89kDa
GeneID?Human?	SLC9A2 NHE2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>solute carrier family 9 member A2(SLC9A2) Homo sapiens This gene encodes a member of the sodium-hydrogen exchanger (NHE) protein family. These proteins are involved in sodium-ion transport by exchanging intracellular hydrogen ions to external sodium ions and help in the regulation of cell pH and volume. The encoded protein is localized to the apical membrane and is involved in apical absorption of sodium. [provided by RefSeq, Jun 2016],</p>